

# HAND HELD PULSEOXIMETER

HS-20A



**Self-induced extremely fast detection**



**Suitable for all Patient groups**



**Real time Data Storage for 10 IDs with 80h storage time**



**Anti-Slip & Anti-Disconnection Design**



## Product specifications.

Display Type	2.3" LCD
Patient groups	Adult, Pediatric and Neonate
Power Supply	Dry Battery (22h) or Lithium battery (16h)
Product weight	180g (device only)
Product Size	72 x 142 x 31 mm
Alarm:	adjustable
<b>Measurement accuracy</b>	
SpO2	70% - 100% ±2 digits 0% - 69%: unspecified
PI	±3bpm
<b>Measurement range</b>	
SpO2	35% - 100%
PI	25 - 250 bpm
Protection Class	IPX2

REF	Description	
DHP01003	Accurate HS-20A	<b>DEVICE</b>
DHP01004	Accurate HS-20A Lithium Version	
DHP02001	Adult soft sensor, 1m	<b>SENSOR</b>
DHP02002	Adult fingerclip sensor, 1m	
DHP02003	Pediatric soft sensor, reusable, 1m	
DHP02004	Pediatric fingerclip sensor, reusable, 1m	
DHP02005	Neonate wrap sensor, reusable, 1m	
DHP02006	Reusable SpO2 extension Cable, 2m	
DHP02007	Adult ear clip, reusable, 1m	